
Rolyan® Sof•Stretch® Coil Extension Splint (Capener):


A997-74 through A997-77

Rolyan® Sof•Stretch® Long Coil Extension Splint:

A997-70 through A997-73


A

B

C

D

Be sure to note which splint you are using and follow the instructions that apply to that splint. Unless otherwise stated, these instructions apply to both of the splints.

INDICATIONS

- For contractures up to 45° that feel elastic when stretched.

Rolyan Sof•Stretch Coil Extension Splint (Capener)

- For tightness in the PIP joint and surrounding soft tissue from flexor tendon tightness, flexion contractures, and boutonniere deformity.

Rolyan Sof•Stretch Long Coil Extension Splint

- For flexor tendon, joint, and soft tissue tightness in the DIP and PIP joints.

CONTRAINDICATIONS

- Not for fixed bony contractures.

PRECAUTIONS

- These splints are to be fitted initially by a healthcare professional who is familiar with the purpose for which they are prescribed. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.
- If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.

- Be careful that the straps are not so tight that they interfere with circulation.

CARE AND CLEANING

Hand wash with mild soap and water. Allow to air dry thoroughly before reapplication.

INSTRUCTIONS FOR USE

A properly applied Rolyan Sof•Stretch Coil Extension Splint (Capener) is shown in A; a Rolyan Sof•Stretch Long Coil Extension Splint is shown in C.

1. The proximal pad has a "P".
2. Attach the self-adhesive strap with the adhesive portion over the distal pad on the dorsal surface of the splint.
3. Slide the splint onto the finger until the proximal pad is over the MCP head on the volar surface. The middle pad should be on the dorsum of the finger just proximal to the PIP joint. The spring coils should be on the lateral aspects of the finger at the level of the PIP.

Rolyan Sof•Stretch Coil Extension Splint (Capener)

4. The distal pad should lie just proximal to the DIP joint, and it should not impair DIP motion (A).
5. For greater comfort, pressure on the soft tissue can be decreased by adding a pad of foam strapping material to the volar surface of the strap, as follows (B):

—Cut a piece of foam that is $\frac{1}{2}$ " (1.3cm) wider than the strap and $\frac{1}{4}$ " (1.3cm) longer than needed to cover the volar strap surface.

—Cut a crosswise slit near each end of the pad.

—Thread the strap through it so the foam is next to the skin.

6. The amount of force provided by the splint can be adjusted by tightening or loosening the strap.

Rolyan Sof•Stretch Long Coil Extension Splint

7. The distal pad should lie over the DIP joint (C).
8. Position the tension strap between the two dorsal pads as shown in D. To apply the tension strap:
 - Position the splint on a flat surface so it rests on the proximal pad.
 - Note that the tension strap has hook in the center of one surface. With the hook up, place the strap under both wires of the splint, between the two dorsal pads.
 - Center the hook between the wires.
 - Wrap the strap ends around the wires and onto the hook, fastening the strap so each end covers half of the length of the hook piece. Do not overlap the strap ends.
9. The tension strap should be on the dorsal surface just distal to the PIP joint (C).
10. The amount of force provided by the splint can be adjusted by tightening or loosening the straps.

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