# Rolyan® Splinting Materials

## Rolyan® Aquaplast® Splinting Materials

Splinting materials are intended to be used for fabrication of custom-molded rigid splints, orthoses, and adaptive equipment.

### PRECAUTIONS FOR FINISHED ORTHOSES

**Splint adjustments are to be made only by a healthcare professional. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.**

**If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**

### PATTERNING

Patterns can be drawn on all splinting materials with a scratch awl before placing the material in the water.

### CUTTING

All splinting materials can be cut at room temperature by using heavy-duty shears, or scoring with a utility knife and bending at the scored mark, then completing the cut with a knife or scissors. If more detailed pattern-making is required, it is recommended that the pattern be transferred to the splinting material to make a rough blank. Then heat the rough blank for one-half the recommended time or until material is soft enough to cut with scissors. This method will reduce rough edges and save time with edge finishing. Light splinting materials can be cut at room temperature using scissors.

### HEATING

**For a chart giving heating information for each material, see "Water Temperature and Time Guide for Rolyan and Rolyan Aquaplast Splinting Materials" on page 8.**

The best method for softening all splinting materials is heated water. Use a heat gun only for spot-heating and making adjustments.

- For hot water heating, fill frying pan or heat pan with water at least 1" (2.5 cm) deep. Heat water to temperature recommended in the chart on page 8 (use thermometer to measure water temperature). For Aquaplast splinting materials, use a heat pan liner in the bottom of the heat pan or frying pan. Hydrocollator® units can be used for large pieces of splinting material. Use a heat pan liner to dip the material into the unit, to ensure that the material does not touch the bottom. Do not leave the splint blank unattended.
- In addition to hot water heating, San-Splint, San-Splint Light, and Aquaplast materials can be heated in a convection oven. Note that dry heating makes the material self-bonding, so no surface preparation is required for bonding. Material that is heated in the oven also becomes hotter than material heated in water; check temperature carefully before applying it to the patient. Heat material

without overlapping or folding on a baking sheet with a non-stick surface in a preheated oven as follows.

—San-Splint: heat at 175°F (80°C) for 3 to 4 minutes.

—San-Splint Light: heat at 175°F (80°C) for 2 to 3 minutes.

—Aquaplast: heat at 195° to 205°F (90–95°C) for up to 10 minutes or until transparent or translucent.

*Note: Overheating splinting materials increases the draping/stretching characteristics; cool slightly before handling it to avoid excess stretching.*

### FORMING

After the pattern is made, the splinting material should be heated in water as indicated in the temperature and time chart. The entire pattern or piece of material should feel soft when removed from the water. When Aquaplast materials become transparent or translucent in a heat pan, the material is soft and ready to form. When the material is soft, remove it from the water and pat it dry with an absorbent towel. (To avoid towel impressions or lint, cover the towel with a pillow case before drying the splinting material.) **Check the temperature of the material before beginning to fabricate the splint on the patient.**

Forming some splints on a patient is easier with the assistance of gravity. When forming a volar forearm splint, for example, place the patient's hand in supination whenever possible, to allow the splinting material to drape into the arches of the hand.

*Technical Tip: To prevent hands from sticking to Original Aquaplast and Original Resilient, use cold cream or mineral oil on hands or dip hands in cold water.*

Use smooth strokes when molding splints.

Original Aquaplast®, Aquaplast-T, Watercolors®, ProDrape®-T, Polyform®, Polyflex II®, and TailorSplint® require minimum shaping with gentle strokes. Repeated pressure is not necessary to retain form while cooling. The controlled stretch of Original Resilient®, Resilient-T, Ezeform®, Synergy®, and San-Splint® allows more aggressive handling. These materials require firm pressure to form the splint. All Aquaplast splinting materials have 100% elastic memory to allow repeated molding for easy and economical splint revisions. Re-soften the material and it will return to its original size and shape.

### COOLING

Cold water, cold packs, cold spray, or cold wraps can be used to speed the cooling process.

*Technical Tip: Exercise band that has been stored in a freezer can also be used to speed the cooling process.*

### ADJUSTMENTS

Change in contour or design of the splint can be accomplished easily by spot-heating the area by pouring hot water over it or using a heat gun. Care should be taken not to overheat the material when using a heat gun—keep the hot air stream of the heat gun moving and do not hold it too close to the material.

### ADHESIVE PROPERTIES

Aquaplast-T, Resilient-T, ProDrape-T, Watercolors, Polyform, Polyform Light, Polyflex II, Polyflex Light and TailorSplint have a protective coating to prevent accidental self-bonding. To bond these materials, remove the coating by scraping the surface or brushing the surface with an acetone material. Spot-heat both surfaces with a heat gun until they become tacky, make sure both surfaces are dry, and then press firmly together. This technique will provide a permanent bond, ideal for adding outriggers and other attachments to the finished splint.

San-Splint and San-Splint Light resist bonding, and using a solvent to prepare the surface for permanent bonding is recommended.

Original Aquaplast, Original Resilient, Ezeform, Ezeform Light, and Synergy do not have a protective coating, and they are self-bonding when both pieces are hot and dry. A solvent is not required—but cleaning the splint surface prior to bonding is recommended.

### EDGE FINISHING

Smooth edges by heating the rough edges of the splinting material with a heat gun or dipping the edge into heated water for 5 to 10 seconds. Smooth the heated edges using finger tips or the palm of the hand—a small amount of hand or body lotion on the finger tips makes smoothing the edges easier.

### CARE AND CLEANING

#### Splinting materials

The shelf life of low-temperature splinting materials is indefinite under the following conditions:

- Store at temperatures between 40° and 90°F (4° and 32°C).
- Store at less than 65% relative humidity.
- Avoid prolonged exposure to light, especially ultraviolet light.
- Avoid exposure to corrosive fumes and ethylene oxide fumes.

#### Fabricated splints

The splint will lose its shape in temperatures over 135°F (57°C). It should be kept away from sources of heat such as ovens, hot water, open flames, and sunny car windows.

The splint can be cleaned with soap and lukewarm water. The straps can be washed with soap and water, but the adhesive portion should not be moistened. The splint and straps should be allowed to dry thoroughly before reapplication of the splint.

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© 1993, Revised 6/94, 7/95, 7/96, 9/96, 2/98, 5/98, 2/01, 9/03
AbilityOne Corporation