PULSTEP SPORECUP® XTR 8896

Einleitung

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PULSTEP SPORECUP® XTR 8896

Introduction

Thank you for purchasing our multifunctional Pedometer! While you are keeping your eyes on the road, this unique product checks how many steps you have taken, how far you have walked or jogged, and how many calories you have burnt. It can also accurately check the pulse rate. It is definitely a valuable asset for joggers!

Key Features

- Displays the walking steps (up to 99.999 steps)
- Displays the walking distance (up to 999.99 miles/km)
- Sports timer : Up to 59'59"99
- Displays the target and actual calories burned in one day (up to 9999,9 kcal)
- Pulse rate checking
- Real time clock (12/24 hr mode)

Standard Mode:

1. Default display is the »step count« & »current time«
2. For pulse rate checking, press [Puls] first, then place your finger on the sensor (on the top right). During measurement, the heart icon will keep on blinking. When it is done, the pulse rate will be shown.
3. Press [Clear] to clear the walking steps, distance and calories burned.
4. Press [Mode] to enter the set up mode.

Time Set, step length and body weight Mode:


1. Press [Set] to enter the time set mode.
2. Press [Up] to toggle between 12 hr or 24 hr, and press [Set] again to confirm.
3. The hour digit will blink, and press [Up] to increase the hour digit, then press [Set] to confirm. The minute digit will blink, press [Up] to confirm. You can press hold [Set] for rapid increase.
4. Then input your step length (from 10 to 150 cm) and body weight (from 10 to 150 kg). Default step length is 60 cm and weight is 60 kg.
5. Unit returns to standard mode automatically 1 minute after the last key pressed.

Count Up Mode:

1. Press [Mode] to enter the sport timer mode.
2. The distance and timer will be displayed on the LCD.
3. Press [Set] to start and press [Set] again to pause the timer.
4. Press [Clear] to reset the timer when stopped.
5. Unit returns to standard mode automatically 1 minute after the last key pressed.

Target burned calories and actual calories

1. Press [Mode], the actual distance taken and calories burned will be displayed.
2. Press [Mode] again, you can set your target calories burned. Press [Set] again, the first digit will blink. Press [Up] to increase it, then press [Set]. The digit will blink, and repeat the same process again until you finished your setup (you can set up to 9999.9 kcal).
3. Unit returns to standard mode automatically 1 minute after the last key pressed.
Note: An individual stride may vary in length. Input the appropriate length by for example measuring the distance of 10 steps and average them by dividing by 10 to get your own step length.


Others

Position the unit correctly

Clip the pedometer on your belt, slacks, or training attire. Maintain a vertical position otherwise incorrect readings may be occurred.

Battery Installation:

1. Open the door by turning anti-clockwise.
2. Insert 1 piece CR2032.
(Note: If not disposed of properly can be harmful. Protect the environment by taking exhausted batteries to authorized disposal stations.)
3. Slide in the rib into the slot next to the battery.
4. Close the battery door by turning clockwise.

Specifications:

Power source: DC 3V (1 x Lithium Battery CR2032)
Battery life: >6 months
Product size: 64 x 49 x 27 mm
Product weight: 36 g

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